

From: bounces+533564.97573500.701046@icpbounce.com on behalf of Separation Science e-Learning <elearning.solutions@sepscience.com>
Sent: Friday, April 15, 2011 12:12 PM
To: Hanchett, James (DPH)
Subject: Today in Separation Science


15 April, 2011

To stop receiving these e-Learning Updates select the
'Unsubscribe' link below the email

John Dolan's
Update

Separation Science e-Learning

HPLC Solutions

 1501815 items in database (showing 15). To help protect your privacy, Outlook prevented automatic download of this picture from the Internet.

ISSUE 66

Enantiomer Separations

Recently, a reader sent a question to our "Ask the Doctor" e-mail asking how to get started with a method that required separation of enantiomers by reversed-phase HPLC. Unfortunately, the problem is a non-starter in this context. Enantiomers cannot be separated by reversed-phase techniques – they require some chirality in the system. John Dolan explains more...

[Click for PDF>>](#)

ISSUE 65

Mobile Phase Lifetimes

A reader recently sent in this e-mail: "I am trying to find information that I can put in an SOP regarding how long a mobile phase can be used before it must be discarded" It would be nice if there were cut-and-dried answer for questions like this. Unfortunately, the answer always seems to start with "it depends." Let's look at some important factors that can influence your decision on how long to use a mobile phase...

[Click for PDF>>](#)

ISSUE 64

Why Acid

Recently, I received an "Ask the Doctor" email from a reader asking why formic acid was specified as an additive for the mobile phase in an HPLC method he was using. Formic or trifluoroacetic acid at 0.1% concentrations are common, especially for LC-MS work. There are a number of reasons for adding an acid at low concentration to the mobile phase. Let's look at two of these: the influence on the column and the sample.

[Click for PDF>>](#)

FEATURED APPLICATIONS

Using PeakView Software with the XIC Manager for Screening and Identification with High Confidence based on High Resolution and Accurate Mass LC-MS/MS
Company: AB SCIEX
[Click to read>>](#)


Fast Analysis of Cefepime and Related Impurities on Poroshell 120 EC-C18
Company: Agilent
[Click to read>>](#)

Advantages of Multidimensional Ion Chromatography for Trace Analysis
Company: Metrohm
[Click to read>>](#)


Analyze More Semivolatiles Samples Per Shift Using Split Injection
Company: Restek
[Click to read>>](#)

Improved Analytical Method for the LC/MS/MS Analysis of Chloramphenicol in Shrimp and Other Marine Food Products using Strata™-X Solid Phase Extraction Cartridges and Kinetex® 2.6 µm Core-Shell Columns
Company: Phenomenex
[Click to read>>](#)

Extraction of Vitamin D Metabolites From Human Plasma Using ISOLUTE® SLE+
Company: Biotage
[Click to read>>](#)


 1501815 items in database (showing 15). To help protect your privacy, Outlook prevented automatic download of this picture from the Internet.

 1501815 items in database (showing 15). To help protect your privacy, Outlook prevented automatic download of this picture from the Internet.

 1501815 items in database (showing 15). To help protect your privacy, Outlook prevented automatic download of this picture from the Internet.


'Chromatography Forum' is now part of Separation Science. Click on the icons to read the latest 'hot topics'...


 Liquid Chromatography

 Gas Chromatography

 Hyphenated Techniques

 Sample Prep

 Data Systems

 **ASK THE DOCTOR**
If you have an analytical question for John Dolan...

[CLICK HERE >>](#)

Follow Separation Science:

 1501815 items in database (showing 15). To help protect your privacy, Outlook prevented automatic download of this picture from the Internet.

Published by Eclipse Business Media Ltd
Frederick House | Princes Court | Beam Heath Way | Nantwich | Cheshire CW5 6PQ | United Kingdom
20 Maxwell Road | #09-17 Maxwell House | Singapore 069113









Copyright © 2011 Eclipse Business Media Ltd. All rights reserved.

To sponsor or advertise in this e-newsletter contact dean.gaines@sepscience.com
For editorial contributions, contact david.hills@sepscience.com

This message was sent from Separation Science e-Learning to james.hanchett@state.ma.us. It was sent from: Eclipse Business Media Ltd, Frederick House, Princes Court, Beam Heath Way, Nantwich, Cheshire CW5 6PQ, United Kingdom. You can modify/update your subscription via the link below.



Unsubscribe

View this message in the iContact Community:  [View message](#)  [Comment on this message](#)  [Receive as RSS](#)
Share this message with others:  [del.icio.us](#)  [Digg](#)  [Reddit](#)  [Facebook](#)  [StumbleUpon](#)